

HI TECH IMPLANT

Born to be different...**the only one !**



Surgical micromotor for implantology with phisiodispenser

Exclusive design and production




Via Vecchia Prov. Luccchese 49/FG
51030 Serravalle Pistoiese • Pistoia • Italy
Tel. +39 0573.994355 • Fax +39 0573.919480
www.bydental.com • info@bydental.com



HI TECH IMPLANT

Born to be different...the only one !

Surgical micromotor for implantology with phisiodispenser



Switching Programs

Switching of the coolant pump

Clockwise/anticlockwise running



Smoothly variable motor control up to the preselected maximum speed.
Straight pressure sensor: no mechanical treatment



The rotor of the pump is protected and it is very easy to attach the irrigation tube

"No-Drops" system: to avoid the problem of the drop to the contra-angle when the flow is stopped



The central unit is over the foot pedal during the carry on.

Available innovative small travel case.



135°C

Maintenance free motor. Guaranteed for more than 500 cycles in autostatic. Male and Female connector to central unit built in aluminium with color code

